LABGIC

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Operations Manual

Liquid Nitrogen Biological Container
液氮生物容器

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目 录 contents

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Table of Contents

Chapter I Overview	1
Chapter II Terms of some products	2
Chapter III product model compilation method	2
Chapter IV Schematic diagram of product structure	3
Chapter V Main product structure and its characteristics	4
Chapter VI Product use notice	5
Chapter VII Product quality assurance scope	8
Chapter VIII Test and calculation methods of product static liquid nitrogen storage period	9
Chapter IX Others	10

Chapter I Overview

Our company has long been devoted to the research, development and system integration of low-temperature biological storage, vacuum insulation and cryogenic treatment, providing complete technical solutions for low-temperature applications, product design and manufacturing of special low-temperature equipment.

Our main products include liquid nitrogen biological containers, large-diameter liquid nitrogen containers, self-pressurized liquid nitrogen containers, cryotherapy devices, liquid nitrogen liquid level alarms, low temperature temperature monitors, liquid nitrogen biological storage systems, nitrogen generators and Various accessories.

My company's product application areas:

Livestock industry: frozen storage of livestock semen and embryos;

Medical industry: frozen storage of human stem cells, plasma, internal organs, etc., and treatment of skin diseases, etc .;

Scientific research insitution: incubators and storage equipment for vaccines, viruses and other microorganisms;

Industrial field: cryogenic treatment of tools, tools, cutting tools, mechanical cold assembly, etc.


Our company has always adhered to the strategy of "taking off with quality and relying on science and technology", taking the lead in establishing a "low-temperature storage R & D center" in the same industry in the country, with a scientific management mechanism and efficient technical team.

Our company pays attention to credibility and pays attention to service. While selling products, it provides users with the most economical and reasonable solutions and the latest industry information at home and abroad.

Chapter II Product terminology

Terminology	Definition
Geometric capacity	The space defined by the geometric dimensions of the inner edge of the inner tank (excluding the neck tube).
Diameter	The inner diameter of the neck tube.
Outer diameter	The outer diameter of the product shell cylinder does not contain ribs.
Height	The shape of the product does not include the maximum length of the cover plug and the lifting barrel along the axial direction of the barrel.
Empty weight	The weight of the product when the air in the inner tank is in normal temperature, with a lid plug, and without the barrel.
Static liquid nitrogen storage period	The product is injected into liquid nitrogen for pre-cooling, after the lifting cylinder is reached to reach thermal equilibrium, then the liquid is filled with liquid nitrogen and placed statically for all days.
Product life	The time interval from when the product leaves the factory until the storage period of liquid nitrogen is below 40% of the factory target.


Chapter III Product model compilation method

Example 1: YDS-30 refer to a liquid nitrogen container with capacity of 30L and diameter of 50mm.

Example 2: YDS-50B-80 refer to a liquid nitrogen container with capacity of 50L and diameter of 80mm, and it is suit for liquid nitrogen transportation.

Chapter IV The product structure diagram

- 1) Handle
- 2) Plug
- 3) Neck Tube
- 4) Adsorbent
- 5) Inner Vessel
- 6) Canister
- 7) Multi-layer Insulation
- 8) Outer Shell
- 9) Nameplate
- 10) Evacuating Nozzle

Chapter V

Product main structure and its characteristics

The container is mainly composed of an outer shell, an inner liner, a neck tube, a multi-layer insulation body, a lifting barrel, etc. (see Figure).

1. The container shell and inner tank are made of aluminum alloy, which has the characteristics of light weight, high strength at low temperature and corrosion resistance.

2. The neck tube is made of glass fiber reinforced plastic, which has the characteristics of high mechanical strength and small thermal conductivity.

3. Multi-layer insulation adopts aluminum foil with excellent reflection performance as the reflection screen, and adopts materials with small thermal conductivity and low outgassing rate as thermal insulation materials to reduce heat radiation.

4. The interlayer between the container shell and the liner is in a high vacuum state to prevent thermal convection of the gas, and an adsorbent with a large adsorption capacity at a low temperature is used to ensure the long-term stability and reliability of the container performance and product life is not less than 5 year.

Chapter VI Product usage instructions

1. Before using the product, the user should open the box to check the appearance quality of the product. If the appearance of the container is found to be defective, within 5 working days after receiving the product, return this information to our company (see the back cover for contact details). Agents and distributors will be subject to the relevant provisions of the agent and distributor agreement.

2. The series of liquid nitrogen biological container products are divided into two types of storage and transportation. The storage type container is mainly used for the static storage of indoor liquid nitrogen and frozen items, and cannot be used for vehicle transportation. In order to meet the transportation conditions, the transport type container has been specially designed for shockproof. It can also be used for transportation under the state of filling with liquid nitrogen, but collision and severe vibration should be avoided.

3. If replenishing the container with liquid nitrogen indoors, please pay attention to the operation of opening doors and windows to prevent severe oxygen deficiency in the operating environment.

4. New containers or containers that have been thawed and rewarmed must be filled with a small amount of liquid nitrogen for pre-cooling before use, and close to thermal equilibrium (not vigorous gasification) before being filled with liquid nitrogen. When filling liquid nitrogen, it is advisable to use a pump or a long tube funnel. The filling tube should be inserted close to the bottom of the container, and there must be a gap in the mouth of the container to allow nitrogen to be discharged. The height of the liquid nitrogen in the container should not exceed the level of the lower end of the neck tube. Avoid touching the human body with liquid nitrogen during operation to avoid frostbite.

5. Within 2 to 3 hours after the container is first placed in a frozen item, it is recommended that the user implement a special person to regularly observe whether there is condensation or frost on the outer surface of the container. If these two conditions occur, it indicates the vacuum degree of the container has deteriorated, which will cause the liquid nitrogen in the container to evaporate in a short time, and the container cannot be used normally. Although the probability of this phenomenon is very small, it may be due to failure to load, unload and transport according to product requirements. However, in order to avoid the loss of frozen items, such observation is necessary.

6. Please note that the container shell has been subjected to the atmospheric pressure of the use and transportation environment. During use or transportation, if a severe collision or heavy pressure occurs, the surface of the container shell will be dented, resulting in product damage or failure.

7. The container is designed for storage and transportation of liquid nitrogen. Liquid oxygen is strictly prohibited.

8. The inner tank of the container is under normal pressure, and it is not allowed to use the method of inflation and pressure or sealing the mouth of the container at will.

9. The neck tube is the channel for liquid nitrogen and frozen items to enter and exit.

10. The container vacuum sealing joint is the key component to maintain the container interlayer vacuum. Once damaged, it will immediately cause the container interlayer vacuum to deteriorate and the product can no longer be used. Therefore, the user should not open the container vacuum sealing joint without authorization.


11. Check the storage amount of liquid nitrogen in the container. You can use the weighing method, or use a liquid level ruler, fine wood or bamboo rod to insert liquid nitrogen according to its frost height (liquid level). Do not insert it with a hollow tube to avoid Liquid nitrogen rushed out of the pipe and hurt people. During use, pay attention to timely replenishing liquid nitrogen. When replenishing liquid nitrogen in the container, the height of its liquid level should not be higher than the lower end plane of the neck tube of the container. When handling and handling frozen items, please take care of them.

12. If the container needs to be cleaned, first wash it with neutral detergent, then rinse it with warm water not higher than 50 °C. After washing, the liner should be dried (preferably dried with hot air, its temperature is not higher than 50 °C) and cooled to normal temperature, then filled with liquid nitrogen.

13. The container should be placed in a cool, ventilated, dry environment; during long-term storage, pay attention to regularly check the liquid surface and the outer surface of the container for condensation and frost.


14. How to open and close the container lid. In order to prevent the accidental opening of the container lid (when not locked) caused by abnormal vibration and collision, the container lid equipped with a buckle of our company adopts a unique self-locking structure design. When opening and closing the container lid, the user should operate as follows:

At the lock of the container lid, insert your fingers (index finger and middle finger) to buckle the edge of the lid to apply and maintain a certain force outward, and then flip up to open the container lid (see Figure). When closing the container lid, a slight “click” sound is heard, indicating that the closure of the container lid is in place, and a self-lock has been formed between the upper lid and the base.

15. The correct extraction method of the canister

Hold the lifting hook of a canister with your hand, lift it vertically about 40 ~ 80 mm, and then rotate the lifting hook 180°, and then vertically lift the entire lifting cylinder out of the container (Figure 3). According to the aforementioned method, all the lifting drums in the container can be extracted one by one.

16. Since this product involves many professional technologies, after failure, the container can only be returned to our company for inspection and maintenance.

Chapter VII Product quality assurance scope

If the liquid nitrogen biological container product produced by the company is used normally, and it is found that the use of the product is indeed affected by the manufacturing of the product within one year from the date of receipt by the user, the following measures shall be taken:

1. When the container is filled with liquid nitrogen for the first time, if a large amount of condensate or frost appears on the surface of the container, our company is responsible for free maintenance or replacement.

2. When the container is tested before the first use, if it is found that the product's static liquid nitrogen storage period is lower than the standard, our company is responsible for free maintenance or replacement.

3. Within one year, when the product appears to be a quality problem that our company manufactures, our company is responsible for checking or repairing the product.

4. If you need to repair or replace, you must ship the container with complete accessories, as well as the product certificate with the same number and a copy of the shopping voucher to the company for inspection.

If one of the following situations occurs, the warranty will become invalid automatically:

1. The user opens the container vacuum sealing joint without authorization, causing damage to the container vacuum sealing joint or damage to the valve core seal.

2. The storage container is used for vehicle transportation when it is filled with liquid nitrogen or with frozen items, resulting in product quality problems.

3. Failure to use the product in accordance with the instructions for use of the product and related regulations requires product damage. If the storage type container is used as a transport type container, the outer shell of the container is depressed, the neck tube is scratched, deformed or leaked by the internal pressure, and the inner tank is corroded.

4. If the product has quality problems, our company is responsible for the maintenance and replacement of this product, and does not deal with liquid nitrogen and preservation matters.

Chapter VIII Test and calculation method of product static liquid nitrogen storage period

1.The static liquid nitrogen storage period of the product is tested by weighing method.

2.The test needs to be conducted at an ambient temperature of $20\text{ }^{\circ}\text{C} \pm 3\text{ }^{\circ}\text{C}$ and normal pressure. At the same time, the test cannot be conducted in a directly ventilated site.

3.The test should be carried out after the liquid nitrogen filling rate of the container is not less than 50% (the filling rate of products under 3L is 100%). Measure the average daily evaporation of product q_m .

4.The maximum weight of the weighing instrument (table scale or electronic digital scale, etc.) used in the weighing method should not be greater than 5 times when the product is filled with liquid nitrogen, and its accuracy should be better than or equal to 1/3000.

5.The product static liquid nitrogen storage period t is calculated by the following formula:

$$t = \frac{m - m_e}{q_m}$$

t —Static liquid nitrogen storage period, d;

q_m —Daily evaporation of the tested product, kg/d;

m —The total weight of the product when it is filled with liquid nitrogen, kg;

m_e —The product empty weight, kg.

Note: The test and calculation of the product's static liquid nitrogen storage period are performed in accordance with GB / T 5458-1997 standard. When the test is not performed according to the above-mentioned ambient temperature and atmospheric pressure, or the minimum division value of the measuring instrument is different from the standard requirements, the calculated value may deviate from the actual static liquid nitrogen storage period of the product.

Chapter IX Others

This product is associated with the product instruction manual, certificate of conformity, etc. Please pay attention to check when open the package. If you have any questions, please contact the sales department of our company within 5 working days.

Beijing Labgic Technology Co., Ltd.

Add:No.9 Yumin Street, Area B of the Airport Industrial Zone,
Shunyi District, Beijing 101318 China
Toll Free:400-600-4213
Website:www.labgic.com